

Through-the-Lens (TTL) Custom Loupes

High Resolution Loupes



FEATURES AND BENEFITS

- Super light weight – Only 38g
- Four point ergonomic customization— TrueFit™ – no standard angles, widths, or loupe heights, or generic working distances
- Break-proof frame lenses – the latest materials science and industry have to offer
- Laser aligned working distance
- True Rx—no prescription manipulation
- Titanium frame, safety frame, or sport frame
- Waterproof
- Military-grade Anti-Scratch coating, Anti-Fog, and Anti-Reflection
- Lifetime Guarantee on frame (Titanium only), optical seal, bond of telescope to frame lens and cracking of frame lens
- 10" - 26" customized working distance

F-Series

F-Series TTL Loupe Models and Prices

MAGNIFICATION	FIELD OF VIEW	WEIGHT	MODEL NO.	WORKING DISTANCE	PRICE
2.5x	4" - 5"	43g	F-2.5 TTL	10" - 26" Customized	\$1095



Q-Optics Custom Through the Lens (TTL) style Loupes are custom made to align to the users facial features. We take multiple measurements and customize your working distance down to the nearest 1/2" or centimeter for our international customers. With our TrueFit™ Measurement System we are sure to get it right. Through the Lens style loupes come in black with your choice of frame style and frame color.

M-Series

M-Series TTL Loupe Models and Prices

MAGNIFICATION	FIELD OF VIEW	WEIGHT	MODEL NO.	WORKING DISTANCE	PRICE
2.5x	4" - 5"	38g	M-2.5 TTL	10" - 26" Customized	\$995
3.0x	2" - 3"	36g	M-3.0 TTL	10" - 26" Customized	\$1045
3.5x	1" - 2"	36g	M-3.5 TTL	10" - 26" Customized	\$1095