

Radiant LED Headlight Systems



The Radiant™ LED light system provides maximum brilliance of LED lighting, with long-lasting battery powered portability. The Radiant™ system comes with everything you need in a stylish and ergonomic package.

- Brightness – Up to 15,000 Foot Candles!
- No Memory – Charge Anytime Battery Pack
- Lithium Ion Battery Pack – twice the longevity of NiMH
- Precision Intensity Controls – saves battery life
- Low Battery Warning
- FITS ALL LOUPE SYSTEMS
- Slim, Small, and Super Light



Radiant™ LED Systems Include:

Battery Power Pack

System Charger

Anti-Curing Filter

LED Headlight

Mounting Clip
(See Chart A on Pg. 9)

General System Specifications:

HEADLAMP	
Intensity	Up to 15,000 foot candles (depending on headlight)
LED Life	50,000 hours
Filter	Anti-Curing
BATTERY PACK	
Battery Life	5 hrs (HIGH) - 10 hrs (MED) - 17 hrs (LOW)
Battery Type	Lithium Ion (no memory recharge)
Weight	5oz
Size	2.5" x 3.25"
CHARGER	
Charge time	Up to 4 hours
Input	Universal 100v – 240v
Plug	US, EU, UK, and AU

Kits includes: Headlamp, battery pack/controller, universal charger, choice of mount from Chart A (Pg. 9)

AVAILABLE COLORS



Blue on Ice



Purple on Ice



Silver on Graphite



Flip-up
Anti-Curing Filter