

Lighted Instruments

Fiberoptic handle style diagnostic probe. May be autoclaved or chemiclaved.


HDP \$135.00



● 4.7 mm end

Fiberoptic handle style micro diagnostic probe. Allows pinpoint precision and access to narrow and tight areas.

HDP-02 \$145.00



● 1.8 mm end

Fiberoptic microscope diagnostic probe for use with endodontic microscopic procedures.

MDP \$154.00



● 0.8 mm end

Fiberoptic Pritchard elevator. Probe can be used on either end of the Pritchard elevator.

PE-FO \$158.00



Fiberoptic mouth mirror. Fiberoptic can be removed and used as diagnostic probe.

MM-FO \$135.00




Fiberoptic Seldin 23. One piece clip is removable and may be used on either end of the Seldin.

SE23-FO \$158.00



Fiberoptic Minnesota retractor.

MR-FO \$148.00



Lighted High Volume Aspirators

Fiberoptic 76R-STRAIGHT comes with a stainless steel high volume aspirator but will attach to most plastic high volume aspirators.


Opening	Fits Into	Product #	Price
10 x 13.5 mm	HVE Cut-Off Valve	76RSTR-FO	\$135.00



overall length 6.25 inches

Fiberoptic straw tip holder. Fits into any standard high volume cut-off valve.

Opening	Fits Into	Product #	Price
7 mm	HVE Cut-Off Valve	STH-FO	\$135.00
Straws for Straw Tip Holder, bag of 100		STH-S	\$9.00



overall length 3.0 inches