

# Fiber Optic Headlight



Mounting Clip  
(See Chart A on Pg. 9)

Headlight

## Fiber Optic Headlight:

Headlight includes:

- Fiber Optic Headlight with universal quick connection
- Light Guide
- Mounting Clip of your choice (Pg. 9)

**HLFOU1 \$450.00**

Light Guide  
(Call for information to connect to your existing light source.)



# Headbands

## Super Lightweight Neoprene Headband

Features 360 degree swivel mount – quickly adjust placement of light spot. Neoprene headband with Velcro latch is both very soft and very durable. Lightweight: 1.8 oz. (51 grams).

**M-HB \$25.00**



Both Headbands feature a 360° Swivel Mount for precise light placement  
(Headbands do not include headlight)

Headlight not included



Velcro Style Latch

Headlight not included

## Heavy Duty Plastic Headband

Features 360 degree swivel mount – quickly adjust placement of light spot. Heavy duty soft plastic band with lever style latch. Neoprene padded front. Attachable face shield option for maximum protection. Weight: 5.1 oz. (145 grams).

**M-HB-P \$65.00**

## Faceshield Kit

Includes 3 clear plastic face shields and the hardware to mount them to the HBHL360-P.

**FS KIT \$15.00**

**FS (Replacement face shield - 3 pack - no hardware) \$12.00**



Lever Style Latch

# Fiber Optics

## Q-Optics® LS-ARC Light Source

Transformer converts automatically to adjust for various voltages. Includes US power cord, other power cords available upon request.

- Intensity: 4200 fc
  - Color: 5500°K (Daylight)
  - Bulb Life: 500 hrs
  - Dimmer: Manual Shutter
  - Single port nose piece
- Dimensions:  
 Length: 4.5 in. (11.43 cm)  
 Width: 3.3 in. (8.38 cm)  
 Height: 3.3 in. (8.38 cm)

**LS-ARC-F (With anti-curing filter) \$750.00**

Replacement bulb for Q-Optics Arc light source. AL-1824, 24W, 60 Volt metal halide arc lamp. Color temperature: 5,500K. Intensity: 4000 fc(43,000 lux)

**AL-1824 \$174.00**

Replacement bulb for LS110-F and LS220 light source. EJA, 150W, 21 Volt Quartz Halogen Dichroic Reflector. Color temperature: 3,400K. Lamp Radiance in Lumens: 4,500.

**EJA (1-9 bulbs) \$32.00 ea.**

**EJA (10 or more bulbs) \$15.00 ea.**



## 110V Fiberoptic Light Source with Anti-Curing Filter

- Intensity: 3200 fc
- Color: 3450°K
- Bulb Life: 40 hrs
- Dimmer: Electrical
- Three port nose piece

**LS110-F (110V) Call for pricing and availability**

## 316 Light Guide

Our 316 light guide is 7' long and constructed with autoclavable silicone tubing. Longer light guides are available at an additional charge. The 316 light guide will fit any QA light source and can be adapted to fit other sources.

**LG316-7 \$195.00**

The male end of the light guide attaches to the LS110-F or LS220 light source. Other light guide ends and adapters are available to fit other light sources.

The female connector end of the light guide attaches to any of the fiberoptic aspirators or instruments shown on pages 11-12.

## Light Guide Holder

The light guide holder helps keep your light guide up and out of the way, keeping it from being damaged. The light guide holder mounts using two wood screws (supplied).

**LG-HOLD \$35.00**

