

# Surgical Aspirators – fits 11mm High Volume Valve



overall length 6.25 inches


Tapered tip helps prevent clogging in the tube.



Surgical aspirator made of 3/16" diameter stainless steel. Tapered tip helps prevent clogging in the tube. Recommended cleaning brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	HVE Cut-Off Valve	<b>15P1A</b>	<b>\$30.00</b>
2.5 mm	HVE Cut-Off Valve	<b>15P2A</b>	<b>\$30.00</b>
3.0 mm	HVE Cut-Off Valve	<b>15P3A</b>	<b>\$30.00</b>

- 1.5 mm
- 2.5 mm
- 3.0 mm
- 4.0 mm
- 4 x 5 mm



overall length 6.25 inches

Right angle surgical aspirator made of 3/16" diameter stainless steel. Tapered tip helps prevent clogging. Excellent for socket and endo suction. Recommended cleaning brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	HVE Cut-Off Valve	<b>15P1-90A</b>	<b>\$30.00</b>
2.5 mm	HVE Cut-Off Valve	<b>15P2-90A</b>	<b>\$30.00</b>
3.0 mm	HVE Cut-Off Valve	<b>15P3-90A</b>	<b>\$30.00</b>



overall length 6.25 inches

Surgical aspirator made of 3/16" diameter stainless steel. Tissue relief slots help prevent grabbing of the tongue, cheek, and soft tissue. Recommended cleaning brush size (Pg. 25): Medium.


Opening	Fits Into	Product #	Price
4 x 5 mm	HVE Cut-Off Valve	<b>16RSBBA</b>	<b>\$30.00</b>



overall length 6.25 inches

Surgical aspirator made of 3/16" diameter stainless steel. Tissue relief slots help prevent grabbing of the tongue, cheek, and soft tissue. Recommended cleaning brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
4.0 mm	HVE Cut-Off Valve	<b>16RSA</b>	<b>\$30.00</b>



overall length 6.25 inches

Right angle surgical aspirator made of 3/16" diameter stainless steel. Tissue relief slots help prevent grabbing of the tongue, cheek, and soft tissue. Excellent for socket and endo suction. Recommended cleaning brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
4.0 mm	HVE Cut-Off Valve	<b>16RS-90A</b>	<b>\$30.00</b>



The fully autoclavable plastic Euro adapter allows you to use U.S. diameter tips (11 mm) in European size valves (16 mm).  
**EURO-A \$5.00**