

Surgical Aspirators – fits Surgical Tubing



15P#T overall length 7.5 inches
15P#AT overall length 8.5 inches

Tapered tip helps prevent clogging in the tube.


- 1.5 mm
- 2.5 mm
- 3.0 mm
- 4.0 mm
- 4 x 5 mm

Surgical aspirators for use with 1/4" surgical tubing. 7-1/2" overall length. Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	15P1T	\$32.00
2.5 mm	Surgical Tubing	15P2T	\$32.00
3.0 mm	Surgical Tubing	15P3T	\$32.00

Same as above but with 8-1/2" overall length.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	15P1AT	\$32.00
2.5 mm	Surgical Tubing	15P2AT	\$32.00
3.0 mm	Surgical Tubing	15P3AT	\$32.00



16RSBBT overall length 7.5 inches
16RSBBAT overall length 8.5 inches

Surgical aspirator for use with 1/4" surgical tubing. 7-1/2" overall length. Excellent for suction in lingual areas.

Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
4 x 5 mm	Surgical Tubing	16RSBBT	\$32.00

Same as above but with 8-1/2" overall length.

Opening	Fits Into	Product #	Price
4 x 5 mm	Surgical Tubing	16RSBBAT	\$32.00



16RST overall length 7.5 inches
16RSAT overall length 8.5 inches

Surgical aspirator for use with 1/4" surgical tubing. 7-1/2" overall length. Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
4.0 mm	Surgical Tubing	16RST	\$32.00

Same as above but with 8-1/2" overall length.


Opening	Fits Into	Product #	Price
4.0 mm	Surgical Tubing	16RSAT	\$32.00



overall length 8.5 inches

Surgical aspirators made of 1/4" stainless steel tubing and are 8 1/2" long. Tapered tips prevent clogging in the tube. Designed for use with 1/4" surgical tubing. Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	46P1AT	\$32.00
2.5 mm	Surgical Tubing	46P2AT	\$32.00
3.0 mm	Surgical Tubing	46P3AT	\$32.00



overall length 7.5 inches

Surgical aspirators made of 1/4" stainless steel tubing and are 7 1/2" long. Tapered tips prevent clogging in the tube. Designed for use with 1/4" surgical tubing. Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	46P1M	\$32.00
2.5 mm	Surgical Tubing	46P2M	\$32.00
3.0 mm	Surgical Tubing	46P3M	\$32.00



overall length 8.25 inches

The No-Clog aspirator features a built-in stylet. Sliding the button forward causes the stylet to quickly eject particles of bone or tissue that may be clogging the tip. When the button is returned to its original position, suction resumes. The adapter end will fit high volume cut-off valves or surgical tubing.

No-Clog surgical aspirators with 1.5, 2.5, or 3.0 mm opening and built-in stylet. Disassembles for easy cleaning. Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	B1T	\$105.00
2.5 mm	Surgical Tubing	B2T	\$105.00
3.0 mm	Surgical Tubing	B3T	\$105.00