

Surgical Aspirators – fits Surgical Tubing



Surgical aspirators designed for use with 1/4" surgical tubing.

overall length 9.0 inches

Opening	Fits Into	Product #	Price
1.7 mm	Surgical Tubing	BE04	\$35.00
2.3 mm	Surgical Tubing	BE05	\$35.00



Surgical aspirators designed for use with 1/4" surgical tubing.

overall length 9.0 inches

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	BE15P1	\$35.00
2.5 mm	Surgical Tubing	BE15P2	\$35.00
3.0 mm	Surgical Tubing	BE15P3	\$35.00
4.0 mm	Surgical Tubing	BE16R	\$35.00




Stainless steel surgical suction tips with built-in frazier style vacuum control. Available in 1.5, 2.5, and 3.0 mm openings.

overall length 9.75 inches

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Tubing	HU1	\$58.00
2.5 mm	Surgical Tubing	HU2	\$58.00
3.0 mm	Surgical Tubing	HU3	\$58.00

Titanium Surgical Aspirators – fits Surgical Tubing



overall length 7.5 inches

Titanium aspirators for use with high volume valve. Made of 100% pure Titanium, not coated or plated. They were specifically designed for use with Titanium implants. Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
2.5 mm	HVE Cut-Off Valve	T-46P2AT	\$140.00
3.0 mm	HVE Cut-Off Valve	T-46P3AT	\$140.00

Surgical Handle and Tips

Surgical aspirator for use with surgical handle. Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Handle	14P1	\$18.00
2.5 mm	Surgical Handle	14P2	\$18.00
3.0 mm	Surgical Handle	14P3	\$18.00




overall length 4.25 inches



overall length 5.4 inches

Surgical aspirator for use with surgical handle. Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Surgical Handle	14P1-90	\$18.00
2.5 mm	Surgical Handle	14P2-90	\$18.00
3.0 mm	Surgical Handle	14P3-90	\$18.00



overall length 4.25 inches


Surgical handle fits 1/4" surgical tubing. Tips for surgical handle shown below. Recommended cleaning brush (Pg. 25): Medium.

Fits Into	Product #	Price
Surgical Tubing	SH	\$49.00

Saliva Ejector Tips

Surgical aspirator for use with the saliva ejector valve. Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Saliva Ejector Valve	14P1-1/4	\$26.00
2.5 mm	Saliva Ejector Valve	14P2-1/4	\$26.00
3.0 mm	Saliva Ejector Valve	14P3-1/4	\$26.00



overall length 4.25 inches




overall length 6.25 inches

Surgical aspirator for use with the saliva ejector valve. Same as the 14P series above, but with an overall length of 6-1/4". Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Saliva Ejector Valve	15P1-1/4	\$26.00
2.5 mm	Saliva Ejector Valve	15P2-1/4	\$26.00
3.0 mm	Saliva Ejector Valve	15P3-1/4	\$26.00
1.5 mm	Saliva Ejector Valve	15P1-90-1/4	\$26.00
2.5 mm	Saliva Ejector Valve	15P2-90-1/4	\$26.00
3.0 mm	Saliva Ejector Valve	15P3-90-1/4	\$26.00

Surgical aspirator for use with the saliva ejector valve. Recommended brush size (Pg. 25): Small & Medium.

Opening	Fits Into	Product #	Price
1.5 mm	Saliva Ejector Valve	14P1-90-1/4	\$26.00
2.5 mm	Saliva Ejector Valve	14P2-90-1/4	\$26.00
3.0 mm	Saliva Ejector Valve	14P3-90-1/4	\$26.00



overall length 4.25 inches