

High Volume Aspirators - Disposable



overall length 6.0 inches

Screen-Vac disposable aspirators feature a built in screen and tissue relief slots to reduce grabbing of the tongue, cheek, and soft tissue.
SV (Box of 100) \$29.95

Disposable Soft-Tip Plastic Screens help prevent grabbing of the tongue, cheek, and soft tissue. Will fit most metal or plastic high volume aspirators.

STS (Bag of 100) \$14.95



Straw Tip Holder features vacuum relief hole on the side allowing for adjustment of the vacuum by covering/uncovering the hole. Straw not included. Recommended brush size: Large.

Opening	Fits Into	Product #	Price
7 mm	HVE Cut-Off Valve	STH	\$67.00



overall length 3.0 inches

Straws for Straw Tip Holder.
STH-S (Box of 200) \$8.50

