

High Volume Aspirators

High volume aspirators feature a flattened surface to aid in retraction, and tissue relief slots to prevent grabbing of tongue, cheek, and soft tissue.

Opening	Fits Into	Product #	Price
6 x 4 mm	HVE Cut-Off Valve	46FSA	\$30.00
8 x 5 mm	HVE Cut-Off Valve	56FSS	\$30.00
9.5 x 6.5 mm	HVE Cut-Off Valve	66FSS	\$30.00
11.5 x 7.5 mm	HVE Cut-Off Valve	76FSSS	\$30.00

Same as above but with overall length of 4-1/4".

Opening	Fits Into	Product #	Price
6 x 4 mm	HVE Cut-Off Valve	44FSA	\$30.00
8 x 5 mm	HVE Cut-Off Valve	54FSS	\$30.00
9.5 x 6.5 mm	HVE Cut-Off Valve	64FSS	\$30.00
11.5 x 7.5 mm	HVE Cut-Off Valve	74FSSS	\$30.00

Back Cut high volume aspirators. Same sizes and features as shown at top of page but features a reverse bevel on the front of the aspirator.

Opening	Fits Into	Product #	Price
6 x 4 mm	HVE Cut-Off Valve	46FSABC	\$30.00
8 x 5 mm	HVE Cut-Off Valve	56FSSBC	\$30.00
9.5 x 6.5 mm	HVE Cut-Off Valve	66FSSBC	\$30.00
11.5 x 7.5 mm	HVE Cut-Off Valve	76FSSSBC	\$30.00

Tongue Guard high volume aspirators. Same sizes and features as shown at top of page but features a built in tongue guard on the front.

Opening	Fits Into	Product #	Price
6 x 4 mm	HVE Cut-Off Valve	46FSATG	\$42.00
8 x 5 mm	HVE Cut-Off Valve	56FSSTG	\$42.00
9.5 x 6.5 mm	HVE Cut-Off Valve	66FSSTG	\$42.00
11.5 x 7.5 mm	HVE Cut-Off Valve	76FSSSTG	\$42.00

High volume aspirator features tissue relief slots to prevent grabbing of tongue, cheek, and soft tissue.

Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
5.5 x 6.5 mm	HVE Cut-Off Valve	46RSA	\$30.00
5.5 x 6.5 mm	HVE Cut-Off Valve	44RSA	\$30.00

High volume aspirator with 5 mm opening.

Recommended brush size (Pg. 25): Medium.

Opening	Fits Into	Product #	Price
5 mm	HVE Cut-Off Valve	46RA	\$30.00

High volume aspirator with 7 mm opening.

Recommended brush size (Pg. 25): Large.

Opening	Fits Into	Product #	Price
7 mm	HVE Cut-Off Valve	56R	\$30.00

High volume aspirator with 10 x 13.5 mm opening.

Recommended brush size (Pg. 25): Extra Large.

Opening	Fits Into	Product #	Price
10 x 13.5 mm	HVE Cut-Off Valve	76R-STR	\$24.00

High volume aspirator with 10 x 13.5 mm opening.

Recommended brush size (Pg. 25): Extra Large.

Opening	Fits Into	Product #	Price
10 x 13.5 mm	HVE Cut-Off Valve	76R	\$24.00

Right angle high volume aspirator features flattened surface to aid in retraction. Does not come with blue rubber sleeves.

Recommended brush size (Pg. 25): Large.

Opening	Fits Into	Product #	Price
10 x 13.5 mm	HVE Cut-Off Valve	65F-90	\$30.00



Package of 6 autoclavable blue rubber sleeves for use with the 65F-90 aspirator.
BRS (Package of 6) \$22.00