

Tubing • Tubing Adapters • Tubing Connectors

High Volume



Smooth-Flex HVE Tubing 1/2" I.D.
S-FLEX \$3.85 ft.



PVC HVE Tubing 1/2" I.D.
7 sided internal shape helps
prevent collapsing.
HVE-PVC \$2.40 ft.

Surgical



Autoclavable silicone surgical
tubing. 1/4" I.D. x 3/32" wall.
SST \$3.81 ft.



Latex surgical tubing.
1/4" I.D. x 3/32" wall.
ST \$1.85 ft.



Latex thick wall surgical tubing.
1/4" I.D. x 1/8" wall.
TWST \$2.25 ft.

Saliva Ejector



PVC Saliva ejector tubing.
7/32" I.D. x 1/16" wall.
SET \$2.25 ft.

Irrigation



Autoclavable Tygon tubing.
3/16" I.D. x 1/16" wall.
TYGON \$4.10 ft.



Autoclavable silicone irrigation
tubing. 2.4mm I.D. x 0.8mm wall.
IRT \$2.55 ft.

The surgical tubing holder
provides a way to keep your
surgical tubing and aspirator
up and out of the way while
not in use. Attaches using
two wood screws (supplied).
ST-HOLD \$35.00



The 1/2" to saliva ejector and 1/2" to
surgical tubing wye connectors allow
you to use high volume and a saliva
ejector or surgical tubing from a single
high volume hose.

1/2-SE WYE \$27.00
1/2-ST WYE \$27.00



The 1/4", 5/16", and 1/2" wye connectors
allows you to use two tubing
connections from a single suction line.
1/4" wye fits saliva ejector tubing, 5/16"
wye fits surgical tubing and 1/2" wye fits
HVE tubing.

1/4 WYE \$8.00
5/16 WYE \$8.00
1/2 WYE \$8.00



Plastic straight tubing connectors allows
you to extend a tubing's length by
connecting two pieces together. Available
in 1/2", 5/16", and 1/4".

PTC-1/2 \$8.00
PTC-5/16 \$8.00
PTC-1/4 \$8.00



The stainless steel tubing connector lets
you extend the length of surgical tubing
or connect one type of surgical tubing
to another.

TC \$10.00



The tubing adapter set includes a stainless steel adapter and five
feet of surgical tubing which allows you to connect the surgical
tubing to your existing high volume valve. The tubing adapter set
provides greater wrist flexibility and is lighter weight compared to
a high volume valve. The tubing adapter set is also available with
autoclavable silicone tubing.

TAS (Adapter and 5 ft. latex tubing.) \$31.50
TASS (Adapter and 5 ft. autoclavable silicone tubing.) \$37.00



The Cut-Off Clamp allows you to
effectively cut off the suction from
surgical or tygon tubing by opening or
closing the clamp
COC \$7.00



Tubing adapter used to connect
1/4" surgical tubing to a high
volume cut-off valve.
TA \$24.00




Tubing adapter used to connect
1/4" surgical tubing to a high
volume cut-off valve.
TA-STR \$24.00



Spring clamp allows you to
attach your 1/4" surgical tubing to
your dental unit or chair.
SC \$6.00



The fully autoclavable plastic
Euro adapter allows you to use
U.S. diameter tips (11 mm) in
European size valves (16 mm).
EURO-A \$5.00



The high volume to saliva ejector adapter
allows you to use any aspirator designed
for use in a saliva ejector valve with your
high volume valve. Features a screw-on
removable silicone head.
HVE-SE \$26.00



The surgical tubing to high volume
adapter allows you to use any aspirator
designed for a high volume valve with
surgical tubing. Also provides a quick
disconnect for the aspirators.
ST-HVE \$38.00



In line cut-off valve for use with 1/4" surgical tubing.
One end features a quick disconnect so a small section
of tubing can be quickly removed and sterilized. Offers
the same easy to clean features of our FLIP-SE cut-off
valve. Includes valve and 2 quick disconnects.
ILCOV \$45.00
SE-ST-DISC (Extra disconnect for in line valve) \$8.00